

Report No.: 27129
Date: 12.08.2021
Contact: GB

Task:

Application field: Coal / Power plant / Energy

Material: Charcoal

Feed size: up to 50 mm

Feed quantity: 2400 g (about 10 l)

Material specification(s): moist (residual moisture 5%), medium hard, brittle

Customer requirements(s): Sieve analysis 10 / 20 / 30 mm, ISO 3310/2

Subsequent analysis: granulation test

Solution

Selected Instrument(s): Horizontal Sieve Shaker AS 400 control
Vibratory Sieve Shaker AS 450 basic

Configuration(s) Item nos.:

- 1 x AS 400 control, 100-240 V, 50/60 Hz, incl. test report acc. EN 10204 2.2
- 1 x Clamping device "comfort" for test sieves 400 mm Ø
- 1 x Collecting pan, stainless steel, 400 mm Ø, height 65 mm
- 1 x Test Sieve, ISO 3310/2, 400x65mm, 10.0mm, stainless steel, perforated plate, round hole
- 1 x Test Sieve, ISO 3310/2, 400x65mm, 20.0mm, stainless steel, perforated plate, round hole
- 1 x Test Sieve, ISO 3310/2, 400x65mm, 31.5mm, stainless steel, perforated plate, round hole
- 2 x Intermediate ring stainless steel, Ø 400x65 mm
- 1 x Evaluation software EasySieve, version 5 / USB stick, incl. adapter USB-RS232

Comparison with vibratory sieving:

- 1 x AS 450 basic, 230 V, 50 Hz incl. clamping device "standard"
- 1 x Clamping device "standard", for AS 450 basic

Please note: Other electrical versions of the instrument(s) are available with different item numbers.

Parameter(s): Revolution speed 220 rpm

Time: 10 and 20 min

The application report is based solely on the processing of the available sample material in the indicated amount. No legal claims shall be derived from this test report.

Subject to technical modification and errors.
© Retsch GmbH - www.retsch.com - lab@retsch.com

page 1/2
12.08.2021

Report No.: 27129
Date: 12.08.2021
Contact: GB

Achieved result(s): See measuring protocol EasySieve
A comparison with vibratory sieving using the AS 450 basic creates identical results.

Remark(s): A complete bag (2.4 kg, 10 l) of charcoal is used for test sieving. To allow a proper sample movement intermediate rings are used on sieves with aperture size > 15 mm.

Recommendation: The Horizontal Sieve Shaker AS 400 control is suitable to analyse the sample material under the above mentioned conditions.

Retsch Sieve Analysis


Company:	[REDACTED]		
User:	GBE		
Result file:	C:\EasySieve\SIEVEDAT\27067_Charcoal\27129_AS400c_10min_003.rdf		
Task file:	27067_Charcoal.afg	GKL file:	27067_Charcoal.gkl
Date:	12.08.2021 12:51	Sieving duration [mm:ss]:	10:00
Sieve shaker:	AS 400 control	No. of revolutions [1/min]:	220.0
Manual input:	Yes		

Material: Charcoal

Sample quantity [g]:	2397.40	Backweighing [g]:	2397.40
Sieve loss [%]:	-0.00		

Sampling:	1 package 2.5 kg	Sample preparation:	No sample preparation
Sieving method:	horizontal AS 400	Test sieve size:	400 x 65 mm
Test sieves according to standard:	DIN ISO 3310-2	Sieving aids:	No sieving aids
Comment:	Intermediate rings on sieves > 15 mm		

Size class	[mm]	p3 [%]	Q3 [%]	Δm [g]	$\Sigma \Delta m$ [g]
<	10.000	12.9	12.9	308.40	308.40
10.000 -	20.000	19.4	32.3	465.50	773.90
20.000 -	31.500	34.4	66.7	824.20	1598.10
> 31.500		33.3	100.0	799.30	2397.40



Characteristics:

Q3 [%]	x [mm]
10.0	7.774
50.0	25.927
90.0	42.526
x [mm]	Q3[%]
0.400	0.5
Span:	1.340
U:	3.766
Sv [mm ² /mm ³]:	0.36
Sm [cm ² /g]:	3.6
XSt [mm]:	16.531
AFS-Nr.:	---
CV:	54.145
MA [mm]:	25.927
d'[mm]:	31.342
n:	1.8
Correlation:	0.99

User: _____

Graph:
...AT\27067_Charcoal\27129_AS400c_10min_003.rdf (27067_Charcoal.afg)

