


EVALUATION SOFTWARE

EasySieve pro /
EasySieve pharma



ATION SOFTWARE
EasySieve pro /
EasySieve pharma

YOUR DIGITAL SIEVING COLLEAGUE

Your work deserves efficiency and reliability – EasySieve pro and its variant EasySieve pharma replace our classic EasySieve software – and the functionality has been significantly expanded for you. Control your sieve analyses in a timesaving and error-free manner!

Every day in the laboratory, you are under time pressure. Incorrect sieving wastes valuable working time and delays subsequent processes. EasySieve pro | pharma was developed precisely for your needs: The software facilitates the preparation, execution, and digital evaluation of your particle size analyses and minimizes sources of error – so you can concentrate on the really important tasks.



EasySieve pro / pharma replaces all previous

EASYSIEVE PRO / EASYSIEVE PHARMA

YOUR ADVANTAGES IN APPLICATION

Automatic recording and evaluation of weighing data

No manual transfer of notes entered by hand – the software documents all steps automatically and transparently.

Intuitive operation and modern design

The logical, intuitive user interface can be customized to your workflows. This allows you to maintain an overview even when under time pressure.

Standard-compliant reports at the touch of a button – report editor included!

Meet all quality requirements and save time on documentation. The reports can be designed to your own ideas and/or your company's policy.

Complex transformations presented simply

Graphical and tabular evaluations help you quickly assess your results.

Complete control and data acquisition

Control your sieving machine effortlessly thanks to the guided sieving concept with automatic transfer of weighing data – for reliable, digital processes.

Automatic scale recognition

Common analytical balances are automatically detected and configured – saving you manual settings.

Comprehensive data export

Your results can be reused flexibly – for reports, archiving, or as a basis for further evaluations (Excel, LIMS, etc.).

EasySieve pharma: advanced functions not only for the pharmaceutical industry

The pharma version offers additional options compared to EasySieve pro. Benefit from optional hierarchical user management, integrated audit trail, and encrypted data storage for maximum data integrity.


EASYSIEVE PRO / EASYSIEVE PHARMA

GUIDED SIEVING – STEP BY STEP TO THE DESIRED RESULT

The guided sieving principle of EasySieve pro and EasySieve pharma is your personal digital laboratory assistant. Whether routine or demanding distribution analyses – the software guides you through every step of the process:


Parameter setup and sieve management

- Configuration of all analysis settings
- Sieve management with all available sieves
- Setup of reference curves




Workspace sieving

- Selection of the method
- Sample identification
- Weighing the sieves and the sample
- Carrying out the sieving
- Reweighing the sieves
- Evaluation




Workspace results

- Display of all results of the last sieving - diagram, table, overview
- Print results
- Data export




Workspace comparison

- | Comparison of all measurement data
- | Up to 32 results can be displayed simultaneously
- | Print and export results



Workspace comparison with trend

- | Trend analysis based on various characteristics
- | Fraction based
- | $Q_3(x)$ based



EASYSIEVE PRO / EASYSIEVE PHARMA

ENHANCE YOUR NEW PURCHASE OR UPGRADE YOUR EXISTING DEVICES – WITH EASYSIEVE PRO | PHARMA

You can control the following RETSCH sieving machines with EasySieve pro | pharma:

- ✓ AS 200 / 300 / 400 / 450 control
- ✓ AS 200 jet
- ✓ AS 200 tap



Greater clarity for your laboratory management:
Create profiles for different sieving machines or
balances in the EasySieve pro | pharma interface and
access them as needed. In addition, you can of
course also evaluate the sieving data of other sieving
machines with EasySieve pro | pharma.

WEIGHING ASSISTANT - FAST AND RELIABLE

Be sure to use RETSCH sieves to benefit from our
newly developed Weighing Assistant, which directly
recommends the optimal sieve load. This helps you
avoid incorrect results due to overloading or
underloading and saves you unnecessary repetitions.
A familiar problem: Different sample materials make
correct loading a challenge. EasySieve pro | pharma
takes this worry away from you – after entering the
sample bulk density, the system shows you the
optimal weighing range. Your sieving becomes
reproducible and reliable.

Of course, test sieves from other manufacturers can
be used as usual.

INDIVIDUAL TOLERANCES WITH WEIGH-IN- TOLERANCE


Want maximum control over your weighing? With
the Weigh-In-Tolerance function, you can determine
for yourself how accurately your sample weighing
should be performed. You specify the desired sample
quantity and define the permissible deviations up
and down individually – simply as a percentage. This
allows you to adjust the tolerances exactly to your
specific workflows and ensures that each sample is
analyzed reproducibly. Unlike the Weighing
Assistant, the Weigh-In-Tolerance is based entirely
on your individual process requirements and enables
a particularly precise loading range.

BACKWEIGHING-TOLERANCE FOR EVEN GREATER RELIABILITY

Define reweighing tolerances and receive automatic notifications if your sample does not meet the specified quality limits.

PROACTIVELY AVOID INCORRECT SIEVING – WITH TREND ANALYSIS

Trend analysis allows you to always keep an eye on your sieve analyses. By comparing several sievings using freely selectable parameters, you can identify deviations at an early stage and ensure your quality standards are maintained – easily and reliably.



PROCESS MONITORING

Stay flexible and react quickly: continuous comparisons of your sieving results make changes in the process immediately visible. This allows you to detect process deviations early on and make adjustments directly to ensure error-free production. Good for you: Trend analysis helps you avoid repetition, saves valuable time, and increases the reliability of your analyses.

TEST EQUIPMENT MONITORING

Clogging of the sieve mesh (blinding) and decreasing sieve tension due to aging or intensive use render test sieves unusable. Rely on objective data instead of estimates: Trend analysis shows you immediately when your sieves no longer function in accordance with standards and should be replaced. Regardless of this, we continue to offer you a qualified, standard-compliant inspection of your sieves with the RETSCH calibration service!

EASYSIEVE PRO /
EASYSIEVE PHARMA

ESSENTIAL FEATURES - NOT ONLY FOR PHARMACEUTICAL PROFESSIONALS



GMP READY

Exclusively for EasySieve pharma

Meet all regulatory requirements with ease – including audit trail and traceable user assignment. Your data is stored in encrypted form and cannot be changed. This allows you to maintain an overview at all times and fully trace all processes and user actions. Your sieve analysis is therefore transparent and secure – for consistent traceability of every analysis.

The system meets the software validation requirements according to

- ✓ ISPE GAMP 5
- ✓ FDA 21 CFR Part 11
- ✓ EC GMP Annex 11

for compliance with the highest standards in quality assurance.

COMPLETE EVALUATION

You can choose from a broad variety of evaluation data. Such as:

- | | |
|---------------------------|-------------------------|
| ✓ Various percentiles | ✓ RRSB |
| ✓ Quantiles | ✓ AFS fineness number |
| ✓ Grain size distribution | ✓ Specific surface area |
| ✓ Span value | and much more! |

ADVANCED DATA TRANSFER

Connect your LIMS and transfer data directly to the system – for even greater efficiency and transparency.

EASYSIEVE PRO / EASYSIEVE PHARMA

TECHNICAL DATA

EASY SIEVE PRO

System requirements	Windows 11 or higher
Computer	IBM compatible PC
CPU	1.6 GHz or more, Multicore (recommended Intel i7 or comparable)
RAM	4 GB minimum (8 GB recommended)
Connections	min. 2 x USB 2.0 (sieving machine, balance)
HDD	256 GB
Display resolution	1920 x 1080
Graphical User Interface, guided sieving	Yes
Languages	German and English
Parameter units changeable	size in mm/μm/mesh/inch mass in g/kg time in s/min
SOPs	Yes, integrated in methods
Individual sample labelling	Yes, for each sample, reference to method during data choice and in overview
Parameter setting per test sieve	Yes
LIMS ready	yes
Teamviewer online support	Yes

ADDITIONAL FEATURES EASYSIEVE PHARMA

Software validation	ISPE GAMP 5, EC GMP Annex 11 conform, AuditTrail acc. to FDA 21 CFR part 11
User management	Yes, free configurable levels, individual rights option for each user
Role and password management	Yes
Audit-trail with E-signature	Yes
Enforce audit trail comments	yes
Generation/ Export of user lists and daily reports	Yes
Restore	Yes

www.retsch.com/easysieve

ORDER DATA

EASYSIEVE® EVALUATION SOFTWARE

ALLOWS TO CONTROL ALL RETSCH SIEVE SHAKERS EXCEPT AS 200 JET PRO & AS 200 JET PHARMA
(SUITABLE FOR AS 200 CONTROL B AND AS 300 CONTROL B FROM
FROM YEAR OF MANUFACTURE 2016)

32.645.0034



EasySieve® pro for Windows® 11 (64 bit), English/German

32.645.0035



EasySieve® pharma, incl. AuditTrail & User Management, for Windows® 11 (64 bit),
English/German