## KOGELTRILMOLEN MM 400

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
[Klik om video te bekijken](#)

### Product Video

### THE MOST VERSATILE ALL-ROUNDER OF BALL MILLS

- | Max. speed 30 Hz
- | Horizontal oscillation causes strong impact effects for effective sample processing
- | Up to 8 mm feed size and 5 µm final fineness
- | 2 grinding stations for jars of min 2 ml and max 50 ml, adapter for 10 x 2 ml single use vials and 4 x 50 ml conical centrifugation tubes
- | Steel jars can be precooled manually in liquid nitrogen
- | Calibrated speed and time, small benchtop model, storable SOPs and cycle programs, 7 different jar materials

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
## KOGELTRILMOLEN MM 400

# REPRODUCIBILITY OF MECHANOCHEMICAL REACTIONS IN THE MIXER MILL MM 400

Reproducibility is a fundamental principle of scientific research and is essential for ensuring the credibility and reliability of scientific findings. The Mixer Mill MM 400 was tested regarding the reproducibility within a mechanochemical reaction, and it could be proven that it provides excellent reproducibility during several repetitions, for both clamping positions, and also between different devices. [1]


Minor variations of the frequency from 30 Hz to 29 Hz or 28 Hz have an influence on the yield of the reaction. It is of fundamental interest that the mixer mill maintains a set value, e.g. 30 Hz, and does not deviate from it. A premise which is fulfilled by the MM 400 which comes with a calibration certificate.

The mechanochemical reaction  $\gamma\text{-Al}_2\text{O}_3 + \text{ZnO} \rightarrow \text{ZnAl}_2\text{O}_4$  was conducted for 30 min using 25 ml grinding jars, 2 x 15 mm grinding balls, 1 g educts, at 28 Hz, 29 Hz and 30 Hz five times in a row. The comparison between left and right clamping station showed highly reproducible results, also the comparison between the 5 trials.

XRD patterns after the mechanochemical reaction  $\gamma\text{-Al}_2\text{O}_3 + \text{ZnO} \rightarrow \text{ZnAl}_2\text{O}_4$ : Left: Grinding at 28 Hz, 29 Hz and 30 Hz, results after 5th reaction. Middle: Comparison left and right grinding station at 28 Hz 5th reaction each. Right: Reaction 1 to 5 at 30 Hz, right grinding station. Results presented by the group of Claudia Weidenthaler. [8]

The experiments were repeated using another MM 400 device to compare the results between the two mills. Again, the excellent reproducibility was verified for the 5 tests conducted at 30 Hz, for both, left and right grinding station.

Almost identical results (weight % of educts and product) and reproducibility are obtained with a different MM 400 device. Results presented by the group of Claudia Weidenthaler. [1]

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Kogeltrilmolens worden vaak ingezet voor de homogenisatie van biologische monsters. Het zogenaamde bead beating with small glass beads is an established method for cell disruption of yeasts, microalgae or bacteria. The sample is only moderately warmed in the process which can be reduced to a minimum by pre-cooling.

The MM 400 allows for efficient cell disruption of up to 240 ml cell suspension for DNA/RNA and protein

extraction. For accurate diagnosis of infections, it is possible to isolate intact bacteria from tissue in 8 x 30 ml bottles or 10 x 5 ml vials by using adapters.

The MM 400 can be operated with a range of adapters for single-use vials with the following capacities:

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
## CO-CRYSTAL SCREENING WITH THE MM 400

Co-crystal screening can be effectively performed in Mixer Mills. In a study [9] using the MM 400, 2 ml steel tubes and the corresponding PTFE adapter were employed to co-crystallize theophylline and benzamide in a 1:1 ratio under the following conditions:

- | 60 min milling time
- | 30 Hz frequency
- | One 6 mm steel ball per tube
- | Four experiments without solvent and four with 20  $\mu$ L ethanol

X-ray powder diffraction patterns of the eight resulting samples (shown in green) align closely with the simulated reference pattern of the target co-crystal. All observed signals correspond to the desired product, with no significant additional signals, indicating successful and reproducible co-crystal formation. The MM 400 with 2 ml steel tubes delivers consistent results, and this compatibility extends to the MM 500 series, which can also accommodate 2 ml steel tubes.

XRD patterns after the co-crystal formation of theophylline and benzamide after 60 min milling time in the MM 400 against a simulated reference. Results presented by experiments of Dominik Al-Sabbagh. [2]

## CHEMISTRY IN THE MILL: TEFLON RECYCLING (PTFE) USING MECHANICAL ENERGY


### MECHANOCHEMICAL RECYCLING OF PTFE (TEFLON)

Mechanochemical reactions can also be carried out particularly efficiently using the MM 400. Recent research shows how PTFE (Teflon) can be broken down in the MM 400 through a reaction with sodium, using mechanical energy. The intense movement of the grinding balls provides the energy needed to break the stable carbon-fluorine bonds – without any additional heat or pressure. In this way, large portions of the material can be converted into sodium fluoride and carbon – a promising approach for future recycling processes.

Image on the right: Dr. Erli Lu and Dr. Dominik Kubicki with the Mixer Mill MM 400, which was used to decompose PFAs. [4]

The process using MM 400 was part of the renowned science program – “Forschung aktuell” on Deutschlandfunk. Give it a listen!

**The radio segment is available only in German.**

## FUNCTIONALIZING BIOMASS FOR PHARMA APPLICATIONS VIA MECHANOCHEMISTRY

Mechanochemistry is transforming how functional biomaterials are made, and cationic cellulose is a prime example. Using a solvent-free process, cotton fibers are combined with a catalytic base and a minimal additive, then milled together with the cationic reagent to activate the reaction using the Mixer Mill MM 400. This solid-state approach eliminates water and bulk solvents, dramatically reducing chemical use and waste compared to conventional methods. After milling, a short aging step completes the reaction, delivering highly charged cellulose fibers with exceptional performance. [3]

Optimal reaction conditions: Cotton fibers were milled in a 50 ml stainless steel jar with 3 x 10 mm balls for 5 min at 25 Hz, then EPTMAC was added, and the mixture was milled for additional 30 min. The subsequent aging of the reaction mixture at 50 °C for 24 h, followed by Soxhlet extraction (48 h) and freeze drying, resulted in the isolation of pure cCF material.

### Why is this exciting for pharma?


These cationic fibers show strong electrostatic binding to viruses, enabling efficient removal of pathogens from water and process streams—critical for sterile manufacturing and clean water applications. Beyond filtration, the material offers potential in drug delivery, antimicrobial surfaces, and bioprocessing aids. The process achieves outstanding sustainability metrics aligning with green chemistry principles and industry goals. It also allows precise control over charge density for tailored performance.

This innovation demonstrates how mechanochemistry can deliver high-value, eco-friendly solutions for pharmaceutical production—combining safety, efficiency, and sustainability in one breakthrough approach.

## KOGELTRILMOLEN MM 400 IN-SITU RAMAN SPECTROSCOPY

*In-situ Raman* spectroscopy is a powerful analytical technique that allows for the monitoring and analysis of materials in their natural or process environment. This method utilizes Raman scattering, a phenomenon where light interacts with molecular vibrations, leading to shifts in the wavelength of the scattered light. These shifts provide a unique spectral fingerprint for the material being analyzed, offering insights into its chemical composition or molecular structure.

The "*in-situ*" aspect refers to the ability to observe and measure these characteristics directly during an ongoing process. This can include observing changes in the presence of various chemical reactions, also the so called mechanochemistry. Mechanochemistry involves the use of impact, shearing, or friction actions to induce chemical changes in solids. This approach is increasingly popular for its ability to bypass the need for solvents, potentially offering a more environmentally friendly and energy-efficient pathway for chemical synthesis. The Raman spectroscopy can provide invaluable insights into the reaction mechanism, phase transformations, reaction kinetics or for optimization of reaction conditions.

The MM 400 is "Raman-ready", allowing easy removal of the bottom plate inlay. The bottom plate has openings for the Raman probe to consistently measure at the bottom of the jars by placing the Raman probe underneath the mill and thus underneath the jars, where particle interaction is most intense, ensuring accurate data. The Retsch PMMA grinding jars, with their transparency and chemical resistance, enhance spectral data without contamination. The plane outer shapes of the jars further enhance the spectroscopic data. These design adjustments streamline the experimental workflow. Researchers can now perform *in-situ* Raman spectroscopy with greater ease and precision, opening new possibilities for in-depth material analysis.

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#### 2 ML TUBES FOR CRYOGENIC GRINDING

Small 2 ml steel tubes are used for cryogenic applications. Up to 20 of these tubes can be clamped into the MM 400 using an adapter. The advantage: they can withstand low temperatures and mechanical stress and do not break like disposable vessels. Ideal for the smallest sample quantities in the cryogenic range.

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## REFERENCES

[1] Reaction scheme and performance of the experiments: Prof. Dr. Claudia Weidenthaler, Research Group Leader Heterogeneous Catalysis Powder Diffraction and Surface Spectroscopy, Max-Planck Institut für Kohlenforschung, Mülheim an der Ruhr.

[2] Reaction scheme and performance of the experiments: Dominik Al-Sabbagh, Chemistry Laboratory Technician, Division 6.3 – Structure Analysis, Federal Institute for Materials Research and Testing (BAM), Berlin.

[3] Tatsiana Nikonovich, Yao Yu, Mikko Korkiakoski, Chengji Yang, Iris Seitz, Daniel Langerreiter, Mauri A. Kostiainen, Eduardo Anaya-Plaza, and Sandra Kaabel; Solid-State Synthesis of Cationic Cellulose Fibers from Low-Processed Cotton for Efficient Virus Capture; ACS Sustainable Chemistry & Engineering 2025 13 (42), DOI: 10.1021/acssuschemeng.5c07884

[4] With permission of Dr Erli Lu, Associate Professor in Mechanochemistry & Sustainable Synthesis School of Chemistry, University of Birmingham

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35 ml

22.041.0004

Mixing beakers of polystyrene, 56 ml, 100 pcs.

#### PMMA, TRANSPARENT JARS FOR MECHANOSYNTHESSES

01.462.0539

10 ml, 10 pieces

02.462.0539

10 ml, 2 pieces

#### ACCESSORIES FOR AERATION AND INCREASED PRESSURE

01.462.0548

Jar set incl. aeration jar 28 ml stainless steel, filter 10 µm, sealings and fittings 2 x 1/

8" for hose diameter 3 mm and 0.65 mm wall thickness (2x3 m hose included)

|             |                                                                                                                                       |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------|
| 22.050.0005 | Conversion kit including 2 connectors for the hoses on grinding arm, bottom plate for leading hoses out of the housing, counterweight |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------|

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|             |                                                                                     |                                                                             |
|-------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|
| 22.486.0005 |    | Opening aid for grinding jars, 2 pcs.                                       |
| 02.706.0351 |    | Adapter for use of 2/4 grinding jars 5 ml (01.462.0550)                     |
| 22.085.0007 |    | Gasket for grinding jar 1.5 ml, 10 pcs.                                     |
| 22.085.0008 |    | Gasket for grinding jar 5 ml, 10 pcs. (for grinding jar 01.462.0231)        |
| 22.111.0001 |                                                                                     | Gasket for grinding jar 5 ml, 10 pcs. (for grinding jar 01.462.0550)        |
| 22.085.0009 |  | Gasket for grinding jar 10 ml, 10 pcs.                                      |
| 22.085.0006 |  | Gasket for grinding jar 25 ml hardened steel and stainless steel, 10 pcs.   |
| 22.085.0003 |  | Gasket for grinding jar 25 ml zirconium oxide and tungsten carbide, 10 pcs. |
| 22.085.0005 |  | Gasket for grinding jar 35 ml stainless steel, 10 pcs.                      |
| 22.085.0004 |  | Gasket for grinding jar 35 ml zirconium oxide, 10 pcs.                      |
| 22.085.0002 |  | Gasket for grinding jar 50 ml stainless steel, 10 pcs.                      |

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22.008.0010

Adapter for 5 reaction vials 5.0 ml, made of PTFE

22.008.0014

Adapter for 10 reaction vials 1.5 and 2.0 ml, made of PTFE or stainless steel

22.008.0005

Adapter for 5 reaction vials 1.5 and 2.0 ml, made of PTFE or stainless steel

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22.749.0004

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22.749.0008

Reaction vials made of stainless steel 316L, 2.0 ml, 10 pcs.  
(for use with adapter 22.008.0014)

## MAALKOGELS

GEHARD STAAL

05.368.0029 5 mm Ø

05.368.0030 7 mm Ø

05.368.0059 10 mm Ø

05.368.0032 12 mm Ø

05.368.0108 15 mm Ø

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22.455.0010 2 mm Ø, 500 g (ong. 110 ml)

22.455.0011 3 mm Ø, 500 g (ong. 120 ml)

22.455.0002 3 mm Ø, 200 stuks (ong. 6 ml)

22.455.0001 4 mm Ø, 200 stuks (ong. 14 ml)

22.455.0003 5 mm Ø, 200 stuks (ong. 25 ml)

05.368.0034 5 mm Ø

05.368.0035 7 mm Ø

05.368.0063 10 mm Ø

05.368.0037 12 mm Ø

05.368.0109 15 mm Ø

05.368.0062

20 mm Ø

05.368.0105

25 mm Ø

#### WOLFRAMCARBIDE

22.455.0006

3 mm Ø, 200 stuks (ong. 6 ml)

22.455.0005

4 mm Ø, 200 stuks (ong. 14 ml)

22.455.0004

5 mm Ø, 200 stuks (ong. 25 ml)

05.368.0038

5 mm Ø

05.368.0039

7 mm Ø

05.368.0071

10 mm Ø

05.368.0041

12 mm Ø

05.368.0110

15 mm Ø

#### AGAAT

05.368.0024

5 mm Ø

05.368.0025

7 mm Ø

05.368.0067

10 mm Ø

05.368.0027

12 mm Ø

## ZIRCONIUM OXIDE

|             |                                                                                     |                                |
|-------------|-------------------------------------------------------------------------------------|--------------------------------|
| 05.368.0089 |    | 2 mm Ø, 0.5 kg (ong. 135 ml)   |
| 05.368.0090 |    | 3 mm Ø, 0.5 kg (ong. 140 ml)   |
| 22.455.0007 |    | 3 mm Ø, 200 stuks (ong. 6 ml)  |
| 22.455.0009 |    | 5 mm Ø, 200 stuks (ong. 25 ml) |
| 05.368.0146 |                                                                                     | 7 mm Ø                         |
| 05.368.0094 |    | 10 mm Ø                        |
| 05.368.0096 |    | 12 mm Ø                        |
| 05.368.0113 |  | 15 mm Ø                        |
| 05.368.0093 |  | 20 mm Ø                        |
| 05.368.0106 |  | 25 mm Ø                        |

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|             |                                                                                     |         |
|-------------|-------------------------------------------------------------------------------------|---------|
| 05.368.0045 |  | 10 mm Ø |
| 05.368.0046 |  | 12 mm Ø |
| 05.368.0114 |  | 15 mm Ø |
| 05.368.0047 |  | 20 mm Ø |

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05.368.0042 5 mm Ø

05.368.0043 7 mm Ø

05.368.0044 9 mm Ø

05.368.0003 12 mm Ø

#### GLASPARELS

22.222.0001 0.10 – 0.25 mm Ø, 500 g (ong. 320 ml)

22.222.0002 0.25 – 0.50 mm Ø, 500 g (ong. 320 ml)

22.222.0003 0.50 – 0.75 mm Ø, 500 g (ong. 320 ml)

22.222.0004 0.75 – 1.00 mm Ø, 500 g (ong. 320 ml)

22.222.0005 1.00 – 1.50 mm Ø, 500 g (ong. 320 ml)