
ZEEFTOESTEL AS 200 CONTROL

Worldwide comparability of results

The microprocessor-controlled measuring and control unit of this model ensures a constant vibration height, allowing for 100% reproducibility of results even among different AS 200 control shakers. One particular characteristic makes this RETSCH product stand out from others: Instead of the vibration height, it is possible to set the sieve acceleration which is independent of the power frequency. Together with the possibility of calibration, this ensures comparable and reproducible sieving results worldwide. Thus, all requirements for the test materials monitoring according to DIN EN ISO 9001 are met.

All sieving parameters – vibration height, time, and interval – are set, displayed and monitored digitally which makes operation of the AS 200 control very convenient and quick. Up to 99 standard operating procedures (SOPs) may be stored for routine analyses.

Through the integrated interface the instrument can be connected to a PC and controlled with the evaluation software EasySieve®. This program enables the user to carry out the whole sieving process and its subsequent documentation with convenience, accuracy and conforming to standards.

[Klik om video te bekijken](#)

Product Video

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- | Sieving with 3-D effect
- | For sieves up to (Ø) 203 mm
- | Suitable for dry and wet sieving
- | Measuring range 20 µm to 25 mm
- | Memory for 99 Standard Operating Procedures (SOPs)
- | Digital setting and control of sieving parameters
- | Sieve acceleration independent of power frequency
- | Patented electromagnetic drive (EP 0642844)
- | Test materials monitoring according to DIN EN ISO 9001

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
RETSCH TEST SIEVES AND ACCESSORIES - ENGINEERED FOR SUPERIOR PERFORMANCE

EASYSIEVE / EASYSIEVE CFR EVALUATIE SOFTWARE


EasySieve, the software for particle size analyses, exceeds manual evaluation in many aspects. The software is able to automatically control the necessary measurement and weighing procedures – from the registration of the weight of the sieve up to the evaluation of the data. It is simple and convenient to use and is also available in an FDA 21 CFR Part 11-conform version.

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

ZEEFTOESTEL AS 200 CONTROL

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EFFICIENT POWDER PROCESSING IN DENTAL LABORATORIES WITH THE AS 200 CONTROL

Modern dental laboratories like the Dental Laboratories Rehra in Germany use the vibratory sieve shaker AS 200 control for precise classification of metal powders from 3D printing to ensure consistently high material quality. Reliable sieving ensures that only optimally sized particles are returned to the production process, increasing process stability. Agglomerated and oversized particles are removed in advance. At the same time, efficient powder processing helps to significantly reduce material losses and lower costs.

The reproducible results of the AS 200 control support consistently high quality of dental components such as crowns and implants. Fine blasting media used for surface treatment can also be precisely classified and freed from dust and contaminants. This enables the targeted recovery and reuse of valuable blasting media, improving process consistency and reducing operating costs. A uniform particle size distribution ensures consistent and reliable blasting results.


Metal lasers and 3-D manufacturing in action: Complex structures are built up layer by layer.

Model view after 3-D manufacturing

Oxide layer on a model part before sandblasting

Sandblasting Process for a Model Part

EFFICIENT SEPARATION OF STONES FROM UP TO 9 SOIL SAMPLES WITH THE AS 200 CONTROL

In soil analysis, sieving with the AS 200 control prevents distortion of analytical results, reliably separates rocks from soil samples, and simultaneously minimizes equipment damage. Users achieve time-saving and precise analytical results, while maximizing the service life of their devices.

Unbeatable Together: The AS 200 control in combination with the BB 50 jaw crusher is the proven solution for processing large agglomerates from soil samples and can handle up to nine samples in one batch.

After

A 120 g soil sample with agglomerates up to 15 mm in diameter was processed using a 200 x 25 mm sieve with 2 mm round holes.

APPLICATION AT COLOGNE CATHEDRAL – OPTIMAL PROTECTION FOR HISTORICAL STONE SURFACES

The Cologne Cathedral workshop has relied on the AS 200 control for many years to precisely determine particle size distribution in mortar. Various particle fractions can be accurately separated and combined. The goal is to develop mortar with optimal structure and color that almost perfectly replicates the original stone. The use of the AS 200 control eliminates manual sorting, making mortar production more reproducible and efficient.

The result: mortar with the physical properties to fill cracks and simultaneously prevent water penetration.

BEFORE - AFTER

HIGHEST FOOD QUALITY WITH THE AS 200 CONTROL

At Lebensgarten GmbH, the quality of organic cereal flakes is ensured through precise incoming goods inspections. Sieving analysis separates the flakes into different particle size fractions. Special importance is given to the dust fraction of the product with a particle size $< 500 \mu\text{m}$. This can impair the sealing of packaging on the one hand and negatively affect the consistency of crunchy products on the other.

The patented electromagnetic drive and the three-dimensional throwing motion of the AS 200 control ensure even distribution and optimal separation. Thanks to flexibly adjustable parameters such as amplitude and hole size, various types of flakes can be precisely sieved – ideal for checking raw material quality regarding dust and fine fractions and ensuring high quality in end products.

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99.200.0027

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[LL:iid.retschedlink_test_sieve_range]

ACCESSORIES FOR CLAMPING DEVICES AS 200, AS 300, AS 400

CLAMPING LIDS

32.481.0022

Clamping lid with large window of Perspex for test sieves 200/203 mm Ø

32.481.0014

Universal clamping lid with small window for test sieves 100/150/200/203 mm Ø

32.481.0015

Universal wet sieving lid with small window for test sieves 100/150/200/203 mm Ø

KLEMSTUKKEN

32.142.0001

Clamping nuts, (2 pieces) for clamping device "standard"

32.737.0001

Quick-clamping elements, (2 pieces) for clamping device "comfort" AS 200/300/400

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BEVESTIGINGSSTAVEN

32.248.0002

Threaded rods, (2 pieces) for clamping device "standard"

32.248.0001

Threaded rods, short, (2 pieces) for clamping of max. 5 test sieves for clamping device "standard"

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32.742.0013

Rods, smooth, long (2 pieces) for universal clamping device "comfort",

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ACCESSORIES FOR TEST SIEVES (PANS, RINGS, LIDS)

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