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The AS 300 model has all the benefits of the AS 200 control but is designed for test sieves with a diameter up to 315 mm, providing a sieve surface which is approximately 2.5 times larger. Therefore, the AS 300 is able to separate up to 6 kg of material in one working run. Repetitive operations are greatly simplified with the possibility to store up to 99 standard operating procedures (SOP). For perfectly reproducible sieving results, the AS 300 control can be programmed with sieve acceleration independent of the power frequency instead of vibration height.

The microprocessor-controlled measuring device monitors and automatically readjusts the vibration height. All sieving parameters are set, displayed and monitored digitally. The AS 300 control can be calibrated, and is thus suitable for test materials monitoring. Like all instruments of the "control" series, the AS 300 has an integrated interface for using the evaluation software EasySieve® to control, set and visualize all parameters, including complete documentation of the sieving process

ACCURACY & EFFICIENCY

- | Sieving with 3-D effect
- | For sieves up to (Ø) 315 mm
- | Suitable for dry and wet sieving
- | Measuring range 20 µm to 40 mm
- | Memory for 99 Standard Operating Procedures (SOPs)
- | Digital setting and control of sieving parameters
- | Sieve acceleration independent of power frequency
- | Reproducible and globally comparable sieving results
- | Short sieving times due to large sieve surface and effective movement
- | Test materials monitoring according to DIN EN ISO 9001

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WET SIEVING WITH VIBRATORY SIEVE SHAKERS

There are many applications for which wet sieving is the best solution, e.g. when the material to be tested is a suspension or when a very fine sample ($< 45 \mu\text{m}$) that tends to agglomerate needs to be sieved. All vibratory sieve shakers from RETSCH can be used for wet sieving. There are special accessories like clamping lids with spray nozzle and collecting pans with outlet available. By placing RETSCH's venting rings between the sieves, air cushions can expand without letting liquid or sample material escape.

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IQ/OQ Documents

We provide IQ/OQ documentation for the „control“ sieve shakers to support IQ/OQ certification by our customers.

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Ultrasonic baths and dryers

Suitable for thorough cleaning of test sieves and for quick, gentle drying of samples and sieves.

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RETSCH TEST SIEVES AND ACCESSORIES - ENGINEERED FOR SUPERIOR PERFORMANCE

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EasySieve, the software for particle size analyses, exceeds manual evaluation in many aspects. The software is able to automatically control the necessary measurement and weighing procedures – from the registration of the weight of the sieve up to the evaluation of the data. It is simple and convenient to use and is also available in an FDA 21 CFR Part 11-conform version.

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RETSCH's vibratory sieve shakers are ideally suited for separation, fractioning and particle size determination of cement clinker, chemicals, coffee, construction materials, fertilizers, fillers, flours, grains, metals powders,

minerals, nuts, plastics, sand, seeds, soils and washing powder.

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
To find the best solution for your sample preparation task, visit our application database.

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VIBRATORY SIEVE SHAKER AS 300 CONTROL

Vibratory Sieve Shaker AS 300 control for test sieves up to 305 mm / 12" Ø
(please order clamping device, test sieves and collecting pan separately)

30.033.0001

AS 300 control

100–240 V, 50/60 Hz

incl. test report acc. to EN
10204 2.2

CLAMPING DEVICES AS 300

max. number of fractions, for test sieves Ø

32.662.0008

Clamping device "standard", 11, 305 mm Ø

32.662.0009

Clamping device "comfort", 11, 305 mm Ø

32.662.0012

Wet sieving clamping device "standard", 11, 305 mm Ø

32.662.0014

Wet sieving clamping device "comfort", 11, 305 mm Ø

SIEVE STACKS AND ACCESSORIES AS 300

60.158.000999

Sieve stack consisting of 7 test sieves (ISO 3310-1), 305 mm Ø, 40 mm height
(630 µm, 1.25 mm, 2.50 mm, 5 mm, 10 mm, 20 mm, 31.5 mm) and collecting pan

60.159.000999

Sieve stack consisting of 7 test sieves (ASTM E11), 305 mm (12") Ø, 40 mm height
(30 mesh, 16 mesh, 8 mesh, 4 mesh, 3/8", 3/4", 1 1/4") and collecting pan

22.243.0002

Rubber disc for sieve plate

99.200.0028

IQ/OQ Documentation for AS 300 control

ACCESSORIES FOR CLAMPING DEVICES AS 200, AS 300, AS 400

CLAMPING LIDS

32.481.0022

Clamping lid with large window of Perspex for test sieves 200/203 mm Ø

32.481.0014		Universal clamping lid with small window for test sieves 100/150/200/203 mm Ø
32.481.0015		Universal wet sieving lid with small window for test sieves 100/150/200/203 mm Ø
02.660.0011		Clamping lid with large window of Perspex for test sieves 305 mm Ø
02.660.0017		Wet sieving lid for test sieves 305 mm Ø

CLAMPING ELEMENTS

32.142.0001		Clamping nuts, (2 pieces) for clamping device "standard"
32.737.0001		Quick-clamping elements, (2 pieces) for clamping device "comfort" AS 200/300/400
05.114.0080		O-ring for quick-clamping element for AS 300 and AS 400, 1 piece

CLAMPING RODS

32.248.0002		Threaded rods, (2 pieces) for clamping device "standard"
32.248.0001		Threaded rods, short, (2 pieces) for clamping of max. 5 test sieves for clamping device "standard"
32.742.0010		Rods, smooth, (2 pieces) for clamping device "comfort" AS 300/400

SIEVING AIDS

32.365.0001		Chain ring for test sieves 200 mm and 203 mm Ø to support horizontal sieving
32.050.0001		Brushes, 3 pieces
32.902.0001		Cubes of polyurethane, 12 x 12 x 12 mm, 10 pieces
32.902.0002		Cubes of polyurethane, 20 x 20 x 20 mm, 10 pieces

32.354.0001  Balls of rubber, 20 mm Ø, 5 pieces

32.354.0002  Balls of agate, 10 mm Ø, 10 pieces

32.354.0004 Balls of steatite, 6 mm Ø, 150 g

TEST SIEVE RACK

32.012.0001  Test Sieve Rack for 10 Test Sieves Ø 200 mm/8", height 50 mm/25 mm

ACCESSORIES FOR TEST SIEVES (PANS, RINGS, LIDS)

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