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The screenshot displays the 'Sieving' software interface, specifically the 'Sample' tab. The interface includes a search bar at the top, followed by tabs for 'Fast sieving', 'Manual sieving', and 'Methods'. Below these are sub-tabs for 'Basic', 'Sample', 'Methodology', 'Data', and 'Device'. The 'Sample' tab is active, showing a 'Define sample information' section with fields for 'Sample material*' (AP42), 'Sample preparation*' (No sample preparation), 'Density' (100 g/cm³), 'Source' (Container), and 'Sampling'. A 'Sample weight' field is set to 1000 g. Below this is a 'Sample tolerance and check' section with a toggle for 'Weight-in tolerance' (On) and a range from -5.00 % to 5.00 %. There is also an 'Off' option for 'Weighing assistant'. A 'Comment' field is at the bottom.

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The screenshot shows the 'Sieving' software interface with the 'Sample' tab selected. The interface includes a search bar, navigation buttons (Back, Forward, Home, etc.), and a sidebar with tabs: Basic, Sample, Methodology, Data, and Device. The 'Sample' tab contains the following settings:

- Sample preparation:** Dividing by (Dropdown), Density (100 g/cm³), Source (Container), Sampling (With Spoon), Sample weight (10.00 g).
- Sample tolerance and check:** On (checked), Weigh-in tolerance (-5.00 % To 5.00 %), Off (unchecked), Weighing assistant.
- Expected overgrain:** Mean size (0.02 mm), On (checked), Expected overgrain (7.50 g), From (-15.00 % To 15.00 %).
- Expected overgrain:** Mean size (0.02 mm), On (checked), Expected overgrain (4.00 g), From (-15.00 % To 15.00 %).
- Comment:** (Text input field).

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SOFTWARE FEATURES

Graphical User Interface, guided sieving	是
Mode of operation	Fast sieving, manual sieving and methods
Languages	German and English
Parameter units changeable	Size: mm/µm/mesh/inch Mass: g/kg Time: s/min Pressure: mbar/Pa/psi Volume: ml/cl/l/cm ³
SOPs	Yes, integrated in methods, more than 40
Individual sample labelling	Yes, for each sample, reference to method during data choice and in overview
Parameter setting per test sieve	是
LIMS ready	是
possible export format	.xls, .xlsx, .xps, .csv, .pdf, .png, .jpg, .txt, .rtf, html, .mht
Teamviewer online support	是
Screenshot function	是
Back-up function	是
Automatic storage of results	Yes, in machine, auto report function, PDF document or to printer

www.retsch.cn/AS200jet-pro-pharma

ORDER DATA

AIR JET SIEVING MACHINES AS 200 JET PRO / JET PHARMA

Air Jet Sieving Machines for test sieves 203 mm (8") Ø, with automatic vacuum regulation and internal balance, incl. lids for test sieves 1" and 2" height and rubber mallet - please order test sieve and MANDATORY vacuum cleaner (Item No. 22.748.0020 or Item No. 22.748.0021) separately!

30.034.0001		AS 200 jet pro	100–240 V, 50/60 Hz	incl. test report acc. to EN 10204 2.2 evaluation software installed (no EasySieve required)
30.034.0002		AS 200 jet pharma	100–240 V, 50/60 Hz	incl. test report acc. to EN 10204 2.2 evaluation software installed incl. AuditTrail & User Management (no EasySieve required)

VACUUM CLEANER AS 200 JET PRO / JET PHARMA

22.748.0020	Industrial vacuum cleaner ATTIX 33 230V 50/60Hz (other electrical versions available) MANDATORY vacuum cleaner for 30.034.0001 & 30.034.0002 NOT SUITABLE for 32.935.0008
32.524.0017	HEPA filter for ATTIX 33
32.524.0018	Cover for HEPA filter ATTIX 33
32.524.0019	Replacement filter bags for ATTIX 33, 5 pcs.

CYCLONE AS 200 JET PRO / JET PHARMA

22.935.0042	High-performance cyclone including holder, collecting receptacle (250 ml and 500 ml sample bottle) and 1.5 m FDA certified hose with clamps (industrial vacuum cleaner 22.748.0020 or 22.748.0021 recommended)
05.707.0208	FDA certified hose, 1.5 m
02.183.0151	Connector for cyclone

ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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TEST SIEVES WITH ELECTROFORMED SHEET

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ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0047 IQ/OQ Documentation for AS 200 jet pro

99.200.0046 IQ/OQ Documentation for AS 200 jet pharma

RECALIBRATION AIR JET SIEVING MACHINES

The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents

SER-97.001.0021 Calibration of sieving machines AS 200 jet pro & AS 200 jet pharma
at Retsch/Germany, incl. test report acc. to EN 10204 2.2