

Task:

Application field: Biology

Material: Pig's ears (dried)

Feed size: 50 x 120 mm (after manually pre-cutting)

Feed quantity: 70 g per ear

Material specification(s): medium-hard, tough

Customer requirement(s): < 0.5 - 1 mm

Subsequent analysis: Residue analysis

Solution:

Selected instrument(s): Heavy-Duty Cutting Mill SM 2000
As the ZM 200 was discontinued we recommend to use the ZM 300 with the same settings now
Ultra Centrifugal Mill ZM 200

Configuration(s): SM 2000: 6-disc rotor, stainless steel;
bottom sieve 6 mm, stainless steel
ZM 200: Push-fit rotor, 12 teeth, stainless steel;
ring sieve trapezoid holes 0.75 mm, stainless steel

Parameter(s): Revolution speed SM 2000: approx. 700 rpm
Revolution speed ZM 200: 18000 rpm

Time: ca. 2 - 3 min. (total grinding time)

Achieved result(s): Predominantly < 0.5 mm

Remark(s): The preparation was done in two steps
1. pre cutting in SM 2000
2. fine grinding in ZM 200

Recommendation: For sample preparation of tough materials, the Heavy-Duty Cutting Mill SM 2000 and the Ultra Centrifugal Mill ZM 200 are suitable under the above mentioned conditions.

The application report is based solely on the processing of the available sample material in the indicated amount. No legal claims shall be derived from this test report.

Subject to technical modification and errors.

© Retsch GmbH - www.retsch.com - lab@retsch.com

Pictures of the sample



Fig. 1: Original sample size 150 x 160 mm



Fig. 2: Pre-cut sample size 50 x 120mm



Fig. 3: After pre-grinding in SM 2000



Fig. 4: After fine grinding in ZM 200

The application report is based solely on the processing of the available sample material in the indicated amount. No legal claims shall be derived from this test report.

Subject to technical modification and errors.

© Retsch GmbH - www.retsch.com - lab@retsch.com